

Sharp Object(s)	Nature of Material Used with Sharps	Disposal Container
Needles, Blades, Slides, Cover Slips, etc.	Biohazardous or Non-biohazardous	Sharps container
Broken glass, Pasteur pipettes	Biohazardous	Sharps container or chemically decontaminate and autoclave then place in broken glass bin
	Non-biohazardous	Broken glass bin or hard walled container (i.e. labeled for glass only-cardboard box, plastic bucket, etc.)
Serological pipette tips	Biohazardous	Place in the sharps container or pipette disposal box after chemically decontaminating or autoclaving
	Non-biohazardous	Broken glass bin ( <b>glass pipette</b> ) or regular trash ( <b>plastic</b> )